




# Novel Coronavirus(SARS-COV-2) Rapid Antigen Test Kit



Reg.No. (NMPA) : 20203400831

NMPA&ISO13485&CE

## Applicable scenarios



TEST RESULTS				SIGNIFICANCE
PCR	AG	IgM Ab	IgG Ab	
+	+	-	-	Patient may be in the "incubation period" of SARS-CoV-2 infection.
+	+	+	-	Patient may be in the early stages of infection, and the body's immune response first produced the antibody IgM, but no IgG was produced or the IgG content did not reach the detection limit of the diagnostic reagent.
-/+	-/+	+	-	Patient may be in the middle stage of SARS-CoV-2 infection.
-	-	+	+	Patient is in the late stage of infection, but the human body has developed some immunity SARS-CoV-2 (the persistent antibody IgG has been produced)
-	-	-	+	Patients may have been infected with SARS-CoV-2 in the past, but the patient has recovered or the virus in the body has been cleared.
+	+	-	+	Patient may be in a recurrent stage of infection.

Product information	
Package	Card:25/box
Reaction time	15Min
Storage Condition	4 ~ 30°C
Shelf Life	12months
Sample type	Nasopharyngeal swab or oropharyngeal swab
Product composition	Biosafety bag, test card, Sampling Tubes, Sterile swab